



INTEGRATED BEAM SCANNING & ANNUAL QA

IBA Dosimetry

full range of integrated and independent solutions for Quality Assurance (QA), calibration procedures, and services that maximize efficiency and patient safety in Radiation Therapy and Medical Imaging



INTEGRATED BEAM SCANNING & ANNUAL QA

Beam Data Quality ® Fundamental for Patient Safety



1

Highest quality beam scanning results in the perfect beam dataset.




2

Detailed automated check of each scan ensures the perfect beam data set. Only with SMARTSCAN!



3

Correct TPS commissioning - the key for accurate planning and dose calculation.



4

Correct plan delivery and treatment safety. Your peace of mind for all your patients!

High-quality beam data and correct commissioning of your Linac and TPS are the basis for correct RT planning, treatment delivery for patient safety. However, beam commissioning remains a major challenge in most clinics:

Common Pain Points with Beam Commissioning






Your Solution: Efficient and High-Quality Beam Commissioning

IBA Dosimetry² 45 Years of ® Water Phantom Innovations

It all started with the world's 1st computerized water phantom by Manfred Wellhöfer

Introduction of the Blue Phantom family

1st scanning software integrated on a QA platform




1st and only Complete Solution for Automated & Guided Beam Commissioning



SMARTSCAN™
Automated & Guided Beam Commissioning



Integrated Software





Integrated software for your comprehensive and efficient beam scanning, commissioning and annual QA.

myQA Accept is the most trusted scanning and beam data analysis software for the commissioning of your TPS and Linac, as well as for annual QA.

- Easy workflow interface for all IBA Dosimetry water phantoms
- Comprehensive scanning and commissioning software
- Beam Scanning automation with myQA SMARTSCANö
- Build your commissioning reference database to compare with new installations and to pair Linacs

1

Automated Queue Generation

- Fast & automatic data acquisition with predefined queues for all major TPS vendors
- Save time with smart sorting algorithm for optimized scan sequences
- Maximize scanning efficiency with advanced queue sorting, prioritizing and multiple queue editing
- Intuitive setup of user-specific queues
- Automated data export for TPS beam modeling




2

Data Acquisition

- Adjustable scanning parameters for optimized measurements
- 1D, 2D and 3D graphical and geometrical visualization of detector position during scanning
- Automatic scanning speed adaptation: Fastest scanning with optimal resolution through scanning speed adaptation

Adaptive Scan Optimization (ASO)


- Adapted scanning speed for the different profile segments
- Optimized for accuracy in the penumbra area
- Fast continuous scanning where fewer data points are sufficient



3

Data Analysis

- Accurate data analysis via standard and customizable protocols
- 1D gamma analysis tolerances
- Overlaying profiles for quantitative comparison
- Library of mathematical smoothing and interpolation functions
- Fast creation and export of data tables (PDD, TMR, OAR, etc.)



4

Data Handling

- Easy data handling with advanced filtering and sorting
- Fast creation and data export to all major TPS
- Copy & paste to other applications, e.g. MS Excel
- Exchange data with other IBA Dosimetry applications
- Convenient reporting

Sc.	Field (Energy, Size)	Rattness	Symmetry	FieldWidth	Penum	Data Handling
	X5 - 3 x 3 cm	100,3 %	100,0 %	3,06 cm	0,28	Focus...
	X5 - 3 x 3 cm	100,5 %	100,0 %	3,17 cm	0,32	Calc Corrections...
	X5 - 3 x 3 cm	100,8 %	100,0 %	3,36 cm	0,35	Smooth...
	X5 - 3 x 3 cm	101,1 %	100,0 %	3,68 cm	0,36	Edit...
	X5 - 3 x 3 cm	100,9 %	100,0 %	3,96 cm	0,38	Move...
	X5 - 4 x 4 cm	101,0 %	100,0 %	4,10 cm	0,31	Make Systematic...
	X5 - 4 x 4 cm	101,6 %	100,0 %	4,22 cm	0,36	Rotate...
	X5 - 4 x 4 cm	102,0 %	100,0 %	4,45 cm	0,35	Mathematics...
	X5 - 4 x 4 cm	102,5 %	100,0 %	4,88 cm	0,39	Realign...
	X5 - 4 x 4 cm	103,0 %	100,0 %	5,27 cm	0,45	Align Coordinates...
	X5 - 6 x 6 cm	103,1 %	100,0 %	6,10 cm	0,36	Calculate TMR/TMR...
	X5 - 6 x 6 cm	104,3 %	100,0 %	6,29 cm	0,42	Calculate Array...
	X5 - 6 x 6 cm	105,3 %	100,0 %	6,64 cm	0,48	Calculate Isodose...
	X5 - 6 x 6 cm	106,0 %	100,0 %	7,25 cm	0,54	Extract Scan...
	X5 - 6 x 6 cm	106,7 %	100,0 %	7,85 cm	0,60	Simulation To Desk...
	X5 - 8 x 8 cm	105,9 %	100,0 %	8,13 cm	0,39	Weight Factors...

The Blue Phantom Family Designed for Accuracy & Efficiency




Blue Phantom² The 3D Water Phantom

- The most trusted 3D water phantom solution for comprehensive beam scanning and RTPS/Linac commissioning and annual QA
- The high-end water phantom has a modular design and can be configured and upgraded to fit any need and budget
- Upgradable with SMARTSCANö software for scanning automation and guidance
- Compatible with all standard Linacs, Halcyonö, and CyberKnife[®]#



Blue Phantom^{COMPACT 2} The 2D Water Phantom

- 2D high-end water phantom, half the size of Blue Phantom² and minimized weight
- Promotes easy transportation with more efficient use, e. g. for annual checks
- Ideal for satellite hospitals and commissioning service providers
- Optimized for Varian[®] Halcyonö / Ethosö commissioning validation and annual QA [details next page]#




Blue Phantom^{Helix 2} For TomoTherapy[®]

- 3D scanning solution for the TomoTherapy[®] / Radixact[®] System#
- Find more information on next page

Dedicated Beam Scanning Solutions

GO



NEW: SMARTSCANö automated and guided beam commissioning and annual QA now for Varian# Halcyonö

myQA HALO™ for Varian# Halcyonö / Ethosö

myQA HALOö is the dedicated, proven package solution for independent commissioning & validation and beam data collection of the Varian Halcyonö as well as for monthly and annual scans. The package consists of cutting-edge beam scanning components that are released, clinically implemented, and trusted by over 4,000 satisfied users worldwide!




- Integrated Quality Assurance for your Halcyonö with the #1 Beam Scanning software: **myQA Accept**
- **myQA SMARTSCAN** now supports automated and guided beam validation for Varian# Halcyonö / Ethosö. Faster workflow through templates and optimized queues.
- **Blue Phantom Compact** tank designed for fast collection of all needed scans, low weight for minimal Halcyonö couch sag
- Unique **Stealth** reference chamber including special holder compensates for the absence of light field
- 1-minute validated leveling solution

IBA Dosimetry QA solutions were used by the University of Pennsylvania medical physics team to validate the new Halcyonö linear accelerator for release and clinical use:

We had a very good experience using our Blue Phantom solution to validate our Halcyon™ Linac. The setup was very straightforward, and the flexibility of the myQA Accept software was essential ... The bulk of our validation work was done in 3 days.

Chris Kennedy PhD, DABR

Medical Physicist, University of Pennsylvania, USA



Blue Phantom Helix for TomoTherapy® / Radixact®#

Blue Phantom Helix is dedicated for full 3D scanning of the TomoTherapy/Radixact System. Based on the proven Blue Phantom², this water phantom enables fast and accurate commissioning and QA work optimized for TomoTherapy.

- Optimized 3D water phantom for faster scanning
- Efficient measurements & analysis with myQA Accept
- Certified 0.1mm high positioning accuracy and outstanding reliability
- Long-term mechanical stability

All product and company names are trademarks or registered trademarks of their respective holders. Use of them does not imply any relationship, sponsorship, or endorsement between IBA or its products and the owners of these trademarks.

IBA Dosimetry Water Phantom Innovations




Designed by experts for high-performance and long-lasting beam scanning reliability and accuracy.

1-Minute Leveling, Visually Verified.


Intuitive and precise 4-point interactive micro leveling.

- Faster and more accurate than automated setup
- Visual check provides confidence in setup accuracy
- Avoids time needed for redundancy checks required by automated setup procedures
- Enables permanent visual check of the leveling and water surface level to detect water evaporation during longer commissioning times

1



2



For leveling, simply adjust the four alignment pins towards the water surface 1 until the water surface adhesion touches the tips 2.


Certified 0.1 mm Accuracy

- Uncompromised accuracy for your RTPS and Linac commissioning
- Only IBA Dosimetry water phantoms are calibrated and certified to guarantee the highest accuracy and reproducibility of ± 0.1 mm
- Have full confidence in your beam data




Continuous, Long-Lasting Accuracy

- The unique magnetostrictive sensor technology provides continuous readouts of the water phantom's absolute position in all three axes (even when not moving)
- Certified detector repositioning accuracy of ± 0.1 mm
- The contactless sensor technology minimizes mechanical wear and ensures long-lasting accuracy



Consistent Accuracy in X and Y Axes

- Small ionization chambers like the IBA Dosimetry CC-04 ensure scanning accuracy independent of the scanning direction, regardless of detector movement and orientation (according to TG-106 report)




High-End Clinical Training Courses at IBA Dosimetry

Become an expert in commissioning and beam scanning and learn how to commission with the highest efficiency and data quality!

- Best-practice courses conducted by real dosimetry experts
- Clinical hands-on with SMARTSCAN, myQA Accept and Blue Phantom²
- Special trainings for small fields, FFF, transmission reference chamber, data processing and more



Accessories for Blue Phantoms



Lift Table – Compact Design

- Lift Table with extended legs can be set up without putting weight on the Linac couch ring
- Water Phantom carriage with manual or electric (telescopic) lifting
- Convenient and fast positioning of the water phantom
- Includes leveling frame for vertical and horizontal micro adjustment (electrical version)



Water Reservoir

- Separate tank trolley on wheels with a polyethylene water reservoir
- Small footprint allows convenient storage and easy maneuvering in narrow mazes
- High performance pump for uni-directional or bi-directional water transport to and from the water phantom
- Electronic pump control for TMR/TPR measurement (option)




TMR Set

- For continuous TMR depth dose curve measurement with real-time display of dose vs. water level
- Online measurement of tissue maximum ratio (TMR) with fixed source detector distance
- High-accuracy contactless sensor technology to accurately measure changing water level

Temperature Sensor

- Water temperature measurement in combination with pressure measurement (built-in the CCU)
- ± 0.3 °C measurement accuracy

Beam Data Verification Audit Service




The Challenge: Errors based on poor commissioning data

Errors in beam data collection are reported to frequently lead to flawed TPS beam models and thus to systematic dose computation errors. These systematic beam data errors are difficult to spot and to eliminate¹.

- Commissioning is a challenging process requiring in-depth experience, especially for small field dosimetry.
- Commissioning is usually done under time pressure, and the job is repetitive and error-prone.
- International recommendations for independent audit of the data by a qualified medical physicist (e.g. AAPM TG-106, AAPM TG-53, ESTRO Booklet 10²).

Dosimetric Issues in Radiation Therapy



RT Deficiencies Identified During On-Site Dosimetry®

S. Kry et al.: IJROBP, Vol 99, 5, 2017 P1094-1100

with Monte Carlo

The Solution: Beam Data Verification Audit

To verify the quality and the accuracy of your commissioning beam data you can now benefit from a new and unique audit service based on highly accurate Monte Carlo simulation and analysis.

This service provides a trustful audit result that is independent, high quality and reliable.

The result:


You gain clarity and the peace of mind that the quality of your new or existing beam data is accurate, or know how to improve your data and scans if needed.




- Simply submit your commissioning data for the audit and gain full insight into the quality of your dataset.
- The report includes comprehensive details of the sources of inconsistencies or errors and provides guidance on how to resolve them.
- The audit is available for beam data of any water phantom and for all standard C-Arm Linacs, Halcyon[®], Ethos[®], CyberKnife[®], TomoTherapy[®]/ Radixact[®].

Typical beam data errors detected with Monte Carlo

Wrong Voltage



Penumbra Discrepancy



I have validated our beam commissioning data using the Monte Carlo modeling within SciMoCa. The comprehensive validation report provided guidance to improve the PDD and MLC leaf calibrations in the TPS. We now have the confidence that our beam data is correct and that our treatment plan calculations are of high quality.


M. Kowatsch, Senior Medical Physicist, LKH Feldkirch, Austria

Watch
Webinar



The quality of our measured beam data is critical to ensure accurate treatment planning and delivery, especially for stereotactic and other precision therapies. Medical physicists are challenged by varying responses from diodes and ion chambers because of the inherent differences between the different detectors. For the absolute confirmation of accuracy of beam data, I believe that Monte Carlo validation is the path forward for verifying small field beam data and has broader implications for verification of data for larger field sizes as well.

Raj K. Mitra, PhD, DABR.
Ochsner Health System² New Orleans, LA, USA



DOLPHIN

Advanced Detector for 3D Patient QA

Unique Design

Unique multi-purpose Transmission Detector for Pre-Treatment Patient QA and Machine QA

High-resolution, high-accuracy detector array


- Detector is perpendicular to the beam
- No angular dependence

Optimized for rotational treatments

- Built-in gantry angle sensor
- Time resolved [4D] segment-by-segment response detection

Cable-free design

- Wireless data connection
- Battery-powered; seconds to change




Dolphin accuracy and efficiency in numbers

1513 ionization chambers in Dolphin detector


40x40 cm² full field size for QA of largest fields in one setup

5 mm detector resolution for high accuracy

1 min setup and readiness for QA measurements



Dolphin real Linac output measurements used for 2D and/or 3D plan QA, and for pre-treatment



Check your dose output and profiles, beam uniformity, constancy, flatness, and more

Pre-Treatment QA






Patient pre-treatment QA at the most advanced level

QA efficiency, confidence, and safety before you treat your patients

- Save time for your equipment setup, measurements, and verifications.
- Use Dolphin for pre-treatment QA and Machine QA

|| In our institution we have clinically implemented Dolphin on our three Linacs. Prior to each patient treatment, an advanced measurement-based quality assurance is performed with Dolphin. This is our basis to achieve the highest QA standards and safe treatments for all our IMRT patients today and in the future.




Watch Dolphin user M. Dierl



Patient QA Workflow Simplicity

➤ Easy setup
and QA measurements

➤ Automatic
instant verification results



- Fast detector setup in just 1 minute
- Simply attach Dolphin to your gantry without cables or calibrations
- Transmit your measured data wirelessly for efficient verification at your QA console

- Instant verification of Dolphin measurements vs. TPS plan and display of the result
- “OK” for a test that is successfully passed
- “Evaluate” for discrepancies exceeding your individual threshold
- In case further analysis is required, a full 3D measurement-based evaluation can be performed




Watch 1 min Dolphin setup video



➤ **Optional full 3D patient dose analysis with COMPASS**

The Physicist and Radiation Oncologists **evaluate** the case using IBA's comprehensive 3D / 4D verification software.



- Proven DVH metrics and TPS evaluation tools support comprehensive evaluation and decision making
- Determine when re-planning is required, or when machine-specific tests are necessary
- Understand the clinical impact of the actual delivery in the patient anatomy

Machine QA



myQA Cockpit gives you full overview of your machine QA status and trend analysis. The browser-based application ensures access anytime, anywhere!

Machine QA efficiency

Dosimetry tests with Dolphin Transmission Detector and myQA

- Save time: Check your dosimetry constancy with the same Dolphin setup and in one session with your patient QA
- Easy comparison of test results and trend analysis
- myQA Machines connects Dolphin into your global QA platform



DOLPHIN® & COMPASS

FOCUS ON EVERY PATIENT

DOLPHIN Key Benefits

– Ready to measure in 60 seconds


Simply attach DOLPHIN to the gantry with one „click-in“ at the same aligned and reproducible position with respect to the beam. Along with the cable-free design, this allows operation to start immediately.

– Full field measurement

The DOLPHIN detector design enables measurements at all beam shapes and field sizes: entire large fields up to 40 × 40 cm² at isocenter plane as well as small fields at high precision. As opposed to other solutions, with DOLPHIN the actual patient plan can be measured, which helps prevent additional phantom plans or cumbersome combinations of multiple measurements.

– 3D dose analysis in the patient's CT

The actual measured dose delivery is reconstructed to the patient's anatomy by the COMPASS software, providing a 3D analysis of the actual dose impact. The TPS-class dose engine embedded in COMPASS enables a truly independent second calculation and comparison with the TPS plan data.





Efficient Patient QA - Hardware
Detectors, Phantoms & Accessories

MatriXX - the **ONE** Detector for your Patient QA and Machine QA.

Compatible with:

- ✓ Conventional & IMRT
- ✓ VMAT/RapidArc®
- ✓ FFF Beams
- ✓ Tomo Therapy®



MatriXX^{Evolution}

- ✓ Proven and effective for conventional, IMRT and rotational cases
- ✓ Suitable for high-dose-rate cases (see specifications)


MatriXX^{FFF}

- ✓ Optimized to support current and future high-dose-rate delivery systems
- ✓ Dedicated for high-dose-rate Flattening-Filter-Free (FFF) as well as for conventional measurements

MatriXX Accessories

3D anatomy-based Patient QA with the Gantry Holder

- ✓ Simple and rigid mounting of MatriXX/Startrack
- ✓ Perpendicular to the beam: Optimal for rotational treatments
- ✓ Ideal for 3D Patient QA




MatriXX mounted to the gantry head

Automated gantry angle correction with the Gantry Angle Sensor




More Sensitivity to Fluence Variation

MatriXX ion chambers increase the sensitive area of the array, and the sensitivity to fluence variations by providing the most complete and accurate picture.*




Actual dose at the measurement plane
Graphics shows a zoomed 6cm x 6cm beam across section

Diode Dector



98% of beam information missed. Point measurements limit predictability*

MatriXX Ion Chamber



Accurate and high-sensitivity measurements

Technical Specifications	MatriXX ^{Evolution}	MatriXX ^{FFF}
Number of chambers and type	1020 air-vented ionization chambers	
Read-out time	20ms without dead time (parallel read-out of all chambers)	
Active area	24.4 x 24.4 cm ²	
Sensor layout	MatriXX in a plane arranged in a 32 x 32 grid	
Pixel spacing (center-to-center)	7.6 mm	
Chamber size / volume	4.5 (Ø) x 5 (h) mm ² / 0.080 cm ³	4.5 (Ø) x 2 (h) mm ² / 0.032 cm ³
Nominal sensitivity	2.0nC/Gy	1.4nC/Gy
FFF compatible	Suitable for FFF beams	Optimized for FFF beams
Charge collection efficiency (at 24Gy/min; 10MV FFF at 100cm SDD)	> 97 % at 1.0 mGy/pulse	> 99 % at 1.0 mGy/pulse
Deviation from linearity	≤ 1 % if the dose is ≥ 0.02 Gy	≤ 1 % if the dose is ≥ 0.15 Gy
Temperature & pressure sensor	Automated k(t,p) correction of the chamber signal	
Weight	10 kg. Easy to carry	
Data communication method	Ethernet connection (via standard network cable)	
Gantry Angle Sensor accuracy	+/- 0.6°	



myQA[®] Daily


easy, efficient, and accurate morning QA

myQA[®] Daily

easy, efficient, accurate by design

myQA Daily is the only solution for fast, easy, and high-quality morning Linac QA.

The largest number of ionization chambers provides more beam data for more accurate beam quality verification.



125 ionization chambers

provide the largest amount of measured beam data of any available daily QA device. This means a more accurate morning QA and trend analysis of dose output, flatness, symmetry, center, field size, and energy.

High-resolution centerline measurements

The 31 ionization chambers for each centerline offer greater beam measurement accuracy, especially in the penumbra regions. This allows a finer analysis of daily beam characteristics as well as earlier detection of suspicious trends.

Energy constancy checks

Dedicated ionization chambers with integrated absorber material automatically verify the photon and electron energy constancy "all with the same beam and detector setup. There's no need to manually add buildup material, to change the setup, or to flip the detector.

Field size flexibility

The detector layout provides the flexibility to perform daily QA tests with standard 20x20 cm² or smaller 10x10 cm² beams.

Light field check

Field size markers permit easy verification of the light field's conformity with the radiation field.

Wireless connectivity

The real-time Wi-Fi data exchange and rechargeable battery allow wireless daily QA setup and measurements. The cable-free design enables a convenient workflow and makes it easy to use at multiple Linacs.



Instant results at your fingertips


The server-based myQA Daily software application is the backbone of seamless morning QA checks. The software easily guides the user through just a few steps.

Easy daily QA test interface:

Simply start the application to have your daily test run automatically ready. For test execution, the menu guides you easily through just a few steps.

The server-based software enables easy and fast integration into an existing IT infrastructure. Furthermore, you can execute tests and review your results on any device.

**PROTECT +
ENHANCE +
SAVE LIVES**



Sensor design

Type of detectors	125 ionization chambers, carbon electrodes
Energy verification	Built-in attenuation material
Inherent buildup	6.0 mm ABS

Performance

Photons	Co-60 to 24 MV
Electrons	4 MeV to 24 MeV
Dose	Unlimited
Dose rate	≥ 1 Gy/min, ≤ 24 Gy/min
Dose/pulse	Max. 0.3 cGy/pulse

Connectivity

myQA Daily detector	Ethernet or Wi-Fi connection to the server; additional LAN port
myQA Daily software	Ethernet or Wi-Fi connectivity to the server via the hospital network

Electrical

Power	Battery and ext. battery charger 9V DC power supply [included]
--------------	----------------------------------------------------------------

Software

myQA Daily software/database	Installed centrally on a server or on a PC, web browser application via network workstation, or tablet PC
-------------------------------------	-----------------------------------------------------------------------------------------------------------

DailyQA_Rev1_1018_E | © IBA 2018 | All rights reserved | Manufacturer: IBA Dosimetry GmbH.
 Technical specifications and product features are subject to change without prior notice.
 The device may not be available in your area. For availability, please contact to your local IBA Sales organization.


DETECTORS

For Relative and Absolute Dosimetry

Ionization Chambers and Diode Detectors

Unique Small Field Dosimetry Solutions

Reference Chamber



Stealth^{CHAMBER}


- New «perturbation free» reference signal chamber
- Avoid returning to the LINAC room frequently and repositioning [compared to standard reference chambers]
- Excellent reproducible reference signal quality, even for SRS/SBRT fields

Detectors

RAZOR^{DIODE DETECTOR}

High-performance diode detector


- Chip size: 0.95 x 0.95 mm; t = 0.4 mm
- Sensitive area: \varnothing 0.6 mm
- For Photon and Electron beams in RT



RAZOR^{CHAMBER}

Compact air ionization chamber


- Cavity volume: 0.01 ccm
- Central electrode material: graphite
- For Photon and Electron beams in RT



RAZOR^{NANO CHAMBER}

Smallest available ionization chamber

- Cavity volume: 0.003 ccm
- Central electrode material: graphite
- For Photon and Electron beams in RT



Detectors for Relative and Absolute Dosimetry


IBA Dosimetry offers a full range of ionization chambers and pSi semiconductor detectors for various 2D and 3D water phantom systems and solid phantoms. All detectors are from our in-house production and have been extensively tested to meet the highest criteria in radiotherapy dosimetry.

Air Ionization Chambers

Compact Chambers

Applications:


All compact chambers are designed for measurements with high reproducibility in air, solid, or water phantoms. They are suitable for relative dosimetry of photon, electron and proton fields in radiotherapy.



CC04

CC04

The CC04 is the conventional ionization chambers for measurements of small fields and of ranges with high dose gradients, e.g. stereotactic fields.



CC08 / CC13

CC08


CC08 is used for customized applications during manufacturing and installation of linear accelerators (e.g. TMBuddelship±).

CC13

CC13 is the standard chamber for clinical use in water phantoms and for output factor measurements.

CC25


CC25 is mainly used for quality assurance in air and for low dose measurements in water phantoms.



CC25

CC13-S

The CC13-S is the replacement for the RK chamber used as the standard ionization chamber for RFA phantoms. Parameters of the CC13-S are similar to the CC13.



CC13-S

Farmer Type Chambers

Applications:

All farmer type chambers are designed for measurements with high reproducibility in air, solid, or water phantoms. They are suitable for absolute dosimetry of photon, electron and proton beams in radiotherapy.




FC65-G

FC65-G and FC65-P

FC65-G and FC65-P are the standard reference chambers for clinical reference dosimetry of high-energy photon and electron beams and scientific applications.



FC65-P



FC23-C

FC23-C

FC23-C yields higher precision in measuring of isodose contours.

Because of their robust plastic construction, FC65-P and FC 23-C can be used for daily routine beam check.


Plane Parallel Chambers

Applications:

All plane parallel chambers are designed for measurements with high reproducibility in air, solid, or water phantoms.

PPC05 and PPC40

PPC05 and PPC40 are suitable for absolute dosimetry of electron, photon and proton beams in radiotherapy.



PPC05



PPC40



NACP

NACP

The NACP is designed according to recommendations of the Nordic Association of Clinical Physicists (NACP), Acta Radiologica Oncology 19,55. The chamber is used for absolute dosimetry of electron beams 2-50 MeV. A thinner front wall minimizes contamination of the beam, allows measurements at shallow depths, and guarantees high accuracy even at low electron energies.

Diode Detectors


Applications:

The IBA Dosimetry diode detectors are designed for depth dose and profile measurements in water and air and for output factor measurements in small photon beams.



RAZOR DIODE
especially recommended for SRS/SBRT fields

The diode detectors from IBA Dosimetry are an excellent choice in relative field analysis as well as output factor measurements. They are based on the 3rd generation of pSi semiconductors. Since their introduction in 1992, the high doped p-type silicon detector chips, specifically designed for radiation therapy applications, have been the natural choice for measurements where high spatial resolution is required. The accuracy and lifetime of the diode detectors is unsurpassed in the field of radiation therapy dosimetry.



myQA HALO™

Integrated Quality Assurance for Varian Halcyon™,1


DOSIMETRY

PROTECT +
ENHANCE +
SAVE LIVES

Blue Phantom^{COMPACT}

The right fit for Halcyon™ commissioning and required scans

- Optimized 2D tank design for fast collection of all needed scans
- Fastest tank setup with micro-leveling frame and alignment pins*
- Minimized couch sag through compact and lightweight tank
- Certified to guarantee highest accuracy & reproducibility of min ± 0.1 mm!
- Ideal for your commissioning, as well as other monthly and annual scans



Stealth^{CHAMBER™}

Small Field scanning accuracy with ease

- The Stealth reference chamber is specially designed for accuracy and excellent reproducibility of field sizes down to 0.5 cm x 0.5 cm
- Special holder attaches the Stealth Chamber to the Blue Phantom Compact*

